



SALJOL Rollator LED-Lampe

## Bedienungsanleitung Rollator LED-Lampe

### 01. Detailbeschreibung



### 02. Spezifikation

Spezifikation:  
Material: ABS  
Gewicht: 56g  
Knopf: ON/OFF  
Batterie: 3\*AAA(included)  
Glühlampe: 15 LED  
LED Lebensdauer: 100,000hrs  
Dimension: 7\*3.7cm

### 03. Warnung

Kindern unter 6 Jahren ist die Nutzung nur unter Aufsicht eines Erwachsenen gestattet.

### 04. Gewährleistung & Kundendienst

Gewährleistungszeitraum: 1 Jahr ab Kaufdatum

Im Falle technischer Probleme kontaktieren Sie bitte Ihren Fachhändler oder wenden Sie sich gerne direkt an SALJOL. Wir helfen Ihnen sehr gerne.



## Instruction Manual SALJOL Walker Light

### 01. Detaildescription



### 02. Specification

Specification:  
Material: ABS  
Weight:56g  
Switch:ON/OFF  
Power:3\*AAA(included)  
Bulbs:15 LED  
LED lifespan: 100,000hrs  
Dimension:7\*3.7cm

### 03. Warning

Children under 6 years old should not use this product without adult supervision.

### 04. Warranty & Customer Service

Warranty time: 1 YEAR  
Please don't hesitate to contact the SALJOL customer service if you have questions or problems with the product. We will always try to help you!